


Tag placement is critical to success!

For maximum retention and the lowest risk of bleeding and infection, place ear tags in either of the two spots. (One piece loop tags only fit in the lower spot.) The goal is about 2 inches from the head of the sheep.


The opposite extreme (installing too far from the skull) will result in more lost tags as the tissue on the tip of the ear is tender and tears easily.

Be careful to avoid the large vein (shown at right).

Sheep



Cattle & Goats
(Boer, Kiko & Nubia)



800-282-6631

www.premier1supplies.com

Washington, IA

© Copyright 2014, Premier1Supplies. All rights reserved.